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Dear Customers, thank you for choosing a KaWe product. Our products are known for their high quality and long life. This KaWe product meets the EC Standards and the Low Voltage Directives.

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**Please read these instructions thoroughly and carefully before attempting to use this product and heed the given instructions.**

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**Use:**  
The MedCharge® 3000 charging station is only to be used by authorized professional personnel. The MedCharge® 3000 charging station is to be used exclusively for charging, 2.5 and/or 3.5V NiCd/KaWe-battery charging handles or alternatively NiMH –rechargeable batteries. The charging station recognizes automatically the different battery types (NiCd or NiMH) and controls the charging process correspondingly.

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Charging cells: Type SUB C NiMH 3.5 V with 3000 mAh; type AA NiMH 2.5 V with 1600 mAh or charging cells type AA NiMH 3.5 V with 700 mAh. Type SUB C NiCd 2.5 V with 2600 mAh, type SUB C NiCd 3.5 V with 1300 mAh and type AA NiCd 2.5 V with 1200 mAh as well.

Batterie tipo Sub C NiMH a 3,5 V da 3000 mAh; tipo AA NiMH a 2,5 V da 1600 mAh o tipo AA NiMH a 3,5 V da 700 mAh. Anche tipo Sub C NiCd a 2,5 V da 2600 mAh; tipo Sub C NiCd a 3,5 V da 1300 mAh e tipo AA NiCd a 2,5 V da 1200 mAh.

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**Purpose/operation:**  
After connection of the plug transformer with the main power source, the charging station is ready to be used (the green LED reads „stand by“).

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**Prior to first use, the charging handles are to be charged for a minimum of 10 hours! Caution, do not insert empty batteries into the unit – danger of explosion!**

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After insertion of the charging handle, the LED responsible for the respective charging slot displays the charge state (red = charging, green = charged). When charging is completed, the LED changes from red to green, signifying that the charging handle is completely charged.

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Upon inserting the charging handle into the charging slot, the charging starts automatically. The LED displays the charging process in red. When the charging process is completed, the charging function changes to a trickle charge, which can be seen when the LED becomes green. When the cells to be charged are completely empty, a higher charging current will be used. Thus, the charging time (LED red) will be shortened. In the trickle charge mode, a lower charging current is used, ensuring a permanently charged state of the cells and protection against the build up of excessive heat.

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If so required, for the different charging handles the receiving sleeve in the charging station is to be adapted.

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After unscrewing the cover, the charging battery can be removed and/or be replaced.

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**Unsuited Use / Contraindication of the Products:**  
Any use other than that described here is not in accordance with the intended use of this product. The manufacturer is not liable for any resulting damages. The user alone bears the risk.

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**CAUTION:**  
• Please use only the accessory parts recommended by KaWe. All other connecting parts may cause the risk of fire, electric shocks or injuries.

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• The included alternating current transformer provides a 9V alternating current and 1000 mA. Please do not use any other transformers that do not meet these demands in order to prevent damage to the charging unit or to the transformer.

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• Please do not change the included plug. In the case that the plug is not compatible with the power socket, please have the fitting power socket installed by a qualified technician. Inexpertly installed connections may cause electric shocks.

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• Please do not use the charging unit if it is damaged. If this is the case, it must be repaired by a technician.

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• Please do not disassemble the charging unit. Incorrect reassembly may cause a fire or electric shock.

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• While charging, the batteries and the charging unit may become warm.

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• Please do not use the charging unit as power source for electrical equipment.

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• Protect the charging unit and the transformer against humidity.

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**Capability characteristics:**  
• The charging process is controlled depending on the charge state of the charging cells.  
• Accelerated charge until the end charging voltage is reached.  
• When the battery is full, automatic change-over to trickle charging.  
• Effective overcharge protection.  
• Automatic start of the charging process when the charging circuit is closed, i.e. when the charging handle has been inserted into the charging unit.  
• Galvanic separation between power supply unit and battery by means of transformer in accordance with EN 60335  
• Input voltage 230V +/-10%, 48-62Hz  
• max. main supply 200 mA, sinusoidal.  
• Efficiency factor 90-93%  
• power factor >0.99  
• Max. charging current up to 250 mA  
• Max. charging voltage up to 3.7 V

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**Further information, maintenance, storage:**  
When used and stored properly, the MedCharge® 3000 charging station will serve you many years.

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**Cleaning:**  
The outside of the charging station can be cleaned with a damp, soft and lint-free cloth. To disinfect the unit, a cloth dampened with alcohol can be used. Dipping the unit into solutions or cleaning agents is not allowed!

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**Guarantee:**  
When used under normal circumstances and with attention to these instructions (the users manual), we guarantee this product for up to 2 years from the date of purchase. Should you need further information or should your unit require repair, please contact your dealer!

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**Please dispose of your rechargeable batteries in accordance with the legal regulations for Ni Cd / NiMH**

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**Old electrical equipment: Please dispose of your electric waste in accordance with the legal requirements. Electric waste may not be disposed of together with domestic waste, and instead must be brought to a collection centre where it will be recycled. By complying with this regulation, you not only help protect the environment, but also the health of your fellow men. This regulation applies to the member states of the EU and other European countries with a separate accumulative system for electric equipment.**

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**Common accessory parts:**  
Further information about the MedCharge® 3000 charging station and the charging handles can be attained in the KaWe-price list or in our Internet-Shop at [www.kawemed.de](http://www.kawemed.de)

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**Technical Data of the Charging Station:**  
Input voltage 9V  
Power input: 9V A max.  
Charging current: accelerated charge approx. 250 mA,  
Continuous charging: trickle charge approx. 30 mA  
Service life: indefinite  
Protection class III  
Dimensions: 130/105/70 mm (length/width/height)  
Weight: 375 gr.  
Admissible ambient temperature +10 up to +40°C  
Admissible storage temperature -20 up to +70°C

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**Technical data of the power supply unit:**  
Fuse: integrated thermo-safety fuse  
Service life: indefinite  
Protection class II  
Dimensions 85/75/58 mm (length/width/height)  
Weight: 400 gr.  
Admissible ambient temperature +10 up to +40°C  
Admissible storage temperature -20 up to +70°C

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**Options:**  
- brackets for the wall fixation  
- charging equipment for charging handles of other manufacturers

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