

Serial connection

Serres Suction bags are easily connected in a serial chain for maximum suction capacity. A disposable Serial tubing with a wide-angle connector is installed between the suction bags. Reusable silicone Vacuum tubing connects canisters to the vacuum source. All Serres Suction bags can be used either in single or serial connection use.



57918



57908

Product no.	Product name	Packing size	Reusable
57918	Serial tubing	25	no
57908	Vacuum tubing, silicone	10	yes

In serial connection you can connect up to 6 canisters into one chain. If a Vacuum shift is used to connect two chains, the maximum collection capacity is 36 litres in one trolley. When using Serres Suction bags in series, the liquid flows into the first bag. When it becomes full, the liquid automatically starts to flow into the next bag through the Serial tube without any user intervention. When the last bag in the chain is full or when the procedure is completed, the bags are plugged and disposed of just like in normal use.